
# Primo Star

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**Operating manual**  
**Bedienungsanleitung**  
**Mode d'emploi**  
**Instrucciones de manejo**  
**Manual de instruções**  
**Инструкция по применению**  
**操作手册**  
**دليل المستخدم**

English

Knowledge of this manual is required for the operation of the instrument. Would you therefore please make yourself familiar with the contents of this manual and pay special attention to hints concerning safe operation of the instrument. The specifications are subject to change; the manual is not covered by an update service. Unless expressly authorized, forwarding and duplication of this document, as well as utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation. All rights reserved in the event of granting of patents or registration of a utility model.

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Number of this manual: M60-2-0011 v  
Date of issue: Version 3, 08/29/2008

**CONTENTS**

	Page
<b>1</b>	<b>Notes on Instrument Safety..... 4</b>
1.1	General safety notes..... 4
1.2	Instrument safety and EMC ..... 5
1.3	Unpacking, transportation, storage..... 5
1.4	Disposal..... 5
1.5	Use..... 6
1.6	Warning and information labels..... 8
1.7	Notes on warranty..... 9
<b>2</b>	<b>Description ..... 10</b>
2.1	System overview ..... 10
2.2	Intended use ..... 12
2.3	Instrument description and main features ..... 12
2.4	Objectives..... 13
<b>3</b>	<b>Start-Up and Operation ..... 14</b>
3.1	Starting up the microscope..... 14
3.1.1	Setting up the microscope ..... 14
3.1.2	Mounting the reflected-light fluorescence illuminator ..... 16
3.1.3	Connecting the battery supply unit..... 19
3.1.4	Switching the microscope on / off ..... 20
3.2	Controls of Primo Star (Full-Köhler and Fixed-Köhler) ..... 22
3.3	Controls of Primo Star iLED (Fixed-Köhler) with reflected-light fluorescence illuminator ..... 24
3.4	Operating the microscope ..... 25
3.4.1	Setting interpupillary distance and viewing height ..... 25
3.4.2	Compensating for ametropia and inserting the eyepiece pointer or eyepiece micrometer .... 25
3.4.3	Adjusting transmitted-light brightfield on the Full-Köhler microscope ..... 26
3.4.4	Adjusting transmitted-light brightfield on the Fixed-Köhler microscope ..... 27
3.4.5	Adjusting transmitted-light phase contrast or transmitted-light darkfield..... 28
3.4.6	Adjusting reflected-light fluorescence ..... 30
3.5	Converting the microscope..... 31
3.5.1	Changing the tube ..... 31
3.5.2	Inserting color filters ..... 31
3.5.3	Replacing the 6V 30W halogen lamp or the LED illumination ..... 32
3.5.4	Changing objectives ..... 32
3.5.5	Installing/removing the mirror..... 33
3.5.6	Installing a camera..... 34
<b>4</b>	<b>Care and Troubleshooting ..... 35</b>
4.1	Instrument care ..... 35
4.2	Troubleshooting ..... 36
4.3	Changing the LED module in the reflected-light fluorescence illuminator ..... 37
<b>5</b>	<b>Appendix..... 39</b>
5.1	Technical data ..... 39

## 1 NOTES ON INSTRUMENT SAFETY

### 1.1 General safety notes

Please read this Operating Manual carefully before starting up the microscope.

If you need supplementary information, contact our Carl Zeiss Service or an authorized agency.

To ensure safe operation and troublefree function of the microscope, strictly observe the precautions and warnings given in this manual.

These are marked herein as follows:

**CAUTION**

This symbol indicates a possible hazard to the user of the instrument.

**CAUTION**

Hot surface!

**CAUTION: LED radiation**

LED class 3B, max. 60 mW, 365 - 625 nm

Do not expose yourself to the beam. Avoid radiation exposure to the skin!

**CAUTION**

This symbol indicates a possible hazard to the instrument or system.

**CAUTION**

Disconnect the plug-in power unit from line power before opening the microscope!

**NOTE**

This symbol refers you to advice that you must observe under all circumstances.

## 1.2 Instrument safety and EMC

The Primo Star and Primo Star iLED microscopes have been designed, produced and tested in compliance with the standards DIN EN 61010-1 (IEC 61010-1) and IEC 61010-2-101 "Safety requirements for electrical equipment for measurement, control and laboratory use".

The Primo Star and Primo Star iLED microscopes meet the requirements of the EC Directive 98/79/EC Annex 1 and carry the  mark.

Radio interference suppression in compliance with EN 55011 Class B

Noise immunity in compliance with DIN EN 61326

The instruments have to be disposed of in compliance with the WEEE Directive 2002/96/EC.

## 1.3 Unpacking, transportation, storage

Please observe the following safety notes for unpacking, transportation and storage of the microscope:

- The microscope is supplied packed to commercial standards in a plastic case with cardboard packaging; use the original packaging only for any transportation.
- Retain the original packaging for longer storage or return to the manufacturer.
- When unpacking the equipment, verify that all parts specified on the delivery note are present.

- Keep transportation and storage temperatures as specified in Technical Data.
- Set up the microscope on a stable worktable with solid and smooth tabletop.
- Do not touch optical surfaces to avoid fingerprints.

Risk of burns due to hot surface on the underside of the microscope during operation and up to 10 minutes after power off.

## 1.4 Disposal

Please observe the following safety notes for the disposal of the microscope:

Defective microscopes should not be disposed of with household waste; dispose of them in compliance with the provisions of the law.

The manufacturer of the device is under the legal obligation to take back defective devices.

Batteries of the battery supply unit should not be disposed of with household waste; dispose of them in compliance with the provisions of the law.

## 1.5 Use

The microscopes including their original accessories must not be used for microscopic techniques other than those described in this Operating Manual.

Please observe the following safety notes when using the microscopes:

The manufacturer cannot assume any liability for other applications, including those of individual modules or single components. This also applies to any service or repair work that is not carried out by authorized service personnel. In case of non-compliance, all warranty claims shall be forfeited.

Opening of the device is only allowed to accordingly instructed specialists or to the Service staff.

The LED module of the reflected-light fluorescence illuminator emits LED light of class 3B. For this reason, never look directly into the light and avoid in any case direct exposure of the skin to the light. When working with the microscope, it is absolutely necessary to use the corresponding protective equipment.

Never look into the light beam - neither with nor without optical instruments, not even if you simply wanted to observe the specimen. In case of non-observance your eyes may be damaged!

Do not operate the devices and their accessories included in the delivery in potentially explosive areas, in the presence of volatile anesthetics or combustible solvents, such as alcohol, benzene or similar chemicals.

Dirt and dust may impair the performance of the devices. They must therefore be protected from such influences to the greatest possible extent and covered with the dust cover when not in use. Always check whether the devices are switched off before you cover them (blue power-on light is off).

The microscopes may only be operated by trained personnel who must be aware of the possible dangers involved with microscopy and the particular application concerned. The microscope may be put into operation only if placed on a stable, solid, smooth and hardly flammable surface.

The microscope is a high-precision instrument that can be impaired in its performance or even be destroyed when handled improperly.

The microscope is equipped with a plug-in power unit allowing line voltages to be used in the range between 100 and 240 V  $\pm 10\%$ , 50/60 Hz, without any need for changing the voltage setting on the instrument.

The plug-in power unit meets the requirements of protection class II (with protective insulation). If its casing is damaged, put the plug-in power unit out of operation. The microscope may be operated only with the supplied plug-in power unit.

If it is noted that protection measures are no longer effective, the instrument must be switched off and safeguarded against inadvertent operation. Please contact a Zeiss service agency or the Carl Zeiss Microscopy Service to repair the instrument.

- Always disconnect the power cable before opening the instrument and changing the lamp or LED.

- Wait for the lamp to cool down before replacing it and do not touch the new bulb to avoid fingerprints.

- The instrument may be opened by accordingly instructed specialists or service staff only.
- The operation of the instrument in explosion-risk environments is not allowed.

When using immersion oil, read in any case the safety data sheet.

Immersion oil irritates the skin. Avoid any contact with skin, eyes and clothing. After skin contact, wash the oil off with plenty of water and soap.

After eye contact, immediately rinse the eye with plenty of water for at least five minutes.

If the irritation persists, consult a medical specialist.

Proper disposal of immersion oil: Take care to ensure that immersion oil does not enter surface water or the sewage system.

The microscope is not equipped with special devices for the protection from substances that are corrosive, potentially infectious, toxic, radioactive, or other substances that might be hazardous to health. Make sure to observe all legal regulations, particularly the relevant national accident prevention regulations when handling such substances.

- Before any transportation of the instrument, switch it off and let it cool down (hot surface on the underside of the instrument).

- Operate the device only on a hard, non-combustible base.

- The plug-in power unit must not get in contact with moisture.

1.6 Warning and information labels


Fig. 1 Warning and information labels on the rear side of the stand



Fig. 2 Warning and information labels on the reflected-light fluorescence illuminator

**Fig. 3** Warning and information labels on the battery supply unit

## 1.7 Notes on warranty

The Primo Star and Primo Star iLED microscopes including their original accessories must not be used for microscopic techniques other than those described in this Operating Manual. The manufacturer cannot assume any liability for other applications.

Please note the following information on warranty for the microscopes:


- The manufacturer guarantees that the microscope is free of material or manufacturing defects when delivered.
- Possible defects must be notified to us immediately and steps be taken to minimize damage.
- If notified of such a defect, the manufacturer is obligated to rectify it at its discretion, either by repairing the instrument or delivering an intact replacement.
- No guarantee is provided for defects caused by natural wear (wearing parts in particular) and improper use.
- The instrument manufacturer shall not be liable for damage caused by faulty operation, negligence or any other tampering with the microscope, particularly the removal or replacement of microscope components, or the use of accessories from other manufacturers.


Unauthorized tampering with the instrument shall lead to a forfeit of all warranty claims.

2 DESCRIPTION

2.1 System overview


English

**Case for transport and storage**  
415500-1805-000  
(not suitable for Primo Star iLED)

For Primo Star microscopes  
and Primo Star LED microscope:  
**Battery supply unit**  
415500-1814-000

**Primo Star microscope stand, Fixed-Köhler, R**  
415500-1100-000  
required:  
**Abbe condenser 0.9/1.25 field 18, Fixed-Köhler**  
415500-1701-000  
**Mechanical stage 75x30, drive right**  
415500-1300-000  
**Specimen holder left, for mechanical stage**  
415500-1302-000

**Primo Star microscope stand, Fixed-Köhler, L**  
415500-1102-000  
required:  
**Abbe condenser 0.9/1.25 field 18, Fixed-Köhler**  
415500-1701-000  
**Mechanical stage 75x30, drive left**  
415500-1301-000  
**Specimen holder left, for mechanical stage**  
415500-1302-000

**Primo Star microscope stand, Full-Köhler, R**  
415500-1101-000  
required:  
**Abbe condenser 0.9/1.25 field 20, Full-Köhler**  
415500-1700-000  
**Mechanical stage 75x30, drive right**  
415500-1300-000  
**Specimen holder left, for mechanical stage**  
415500-1302-000


**Primo Star microscope stand, Full-Köhler, L**  
415500-1103-000  
required:  
**Abbe condenser 0.9/1.25 field 20, Full-Köhler**  
415500-1700-000  
**Mechanical stage 75x30, drive left**  
415500-1301-000  
**Specimen holder left, for mechanical stage**  
415500-1302-000

For Full-Köhler only:  
**Slider Ph 1, for Primo Star**  
415500-1815-000  
**Slider Ph 2, for Primo Star**  
415500-1803-000  
**Slider Ph 3, for Primo Star**  
415500-1817-000

**Darkfield slider**  
0.40-0.65  
415500-1802-000

For Full-Köhler only:  
**Slider Ph 1 / H / Ph 2**  
for Primo Star  
415500-1816-000

**Color filter set**  
blue, green, yellow d=45x1.5  
415500-1804-000

**Objective Plan-Achromat 4x/0.10**  
415500-1600-000  
**Objective Plan-Achromat 10x/0.25**  
415500-1601-000  
**Objective Plan-Achromat 10x/0.25 Ph 1**  
415500-1605-000  
**Objective Plan-Achromat 20x/0.40**  
415500-1606-000  
**Objective Plan-Achromat 20x/0.40 Ph 2**  
415500-1607-000  
**Objective Plan-Achromat 40x/0.65**  
415500-1602-000  
**Objective Plan-Achromat 40x/0.65 Ph 2**  
415500-1603-000  
**Objective Plan-Achromat 100x/1.25 Oil**  
415500-1604-000  
**Objective Plan-Achromat 100x/1.25 Ph 3**  
415500-1608-000

**Transmitted-light illuminator**  
with LED  
415500-1200-000

**Transmitted-light illuminator with HAL**  
415500-1201-000  
**Halogen lamp 6 V 30 W**  
415500-1901-000

## 2.2 Intended use

The Primo Star and Primo Star iLED microscopes are universally applicable light microscopes primarily designed for the examination of cell and tissue cultures as well as sediments in culture flasks, Petri dishes and microplates.

Typical applications:

Examination of blood and tissue samples from the human body, observation of intracellular processes on living cell cultures, cell-cell interactions, motility, growth, potential measurements, detection of medical drugs, microinjection and in vitro fertilization.

The Primo Star iLED microscope combined with the reflected-light fluorescence illuminator offers, for example, excellent opportunities for the detection of tuberculosis agents by applying the reflected-light fluorescence method.

When handling hazardous substances, observe the instructions on intended operation, correct use and statutory safety precautions.

## 2.3 Instrument description and main features

The Primo Star and Primo Star iLED microscopes are transmitted-light microscopes of compact design with a small footprint.

Besides the high-resolution, infinity-corrected objectives and the important microscopy techniques, such as brightfield, darkfield and phase contrast in transmitted light, as well as fluorescence in reflected light (only Primo Star iLED with reflected-light fluorescence illuminator), the microscope is optionally available with a camera port for photo and video documentation.

The major features of the microscope include:

- Modular illumination through 6V 30W halogen lamp, LED illumination or illuminating mirror for transmitted light
- Optional reflected-light fluorescence illuminator (Primo Star iLED)
- Integrated fixture for external power supply and cable (incl. cable with multiple plug and country-specific plug inserts)
- Optional rechargeable battery supply unit for network-independent operation of the microscope or for uninterruptible power supply in case of power failure; serves at the same time as a line filter
- Plastic-coated carrying handle integrated in stand for setting up, demounting and transporting the device
- Blue, light-intensity indicators installed on both sides and well visible from a distance
- Primo Star stand in "Full-Köhler" or "Fixed-Köhler" design
- Primo Star iLED stand in "Fixed-Köhler" design
- Convenient coaxial coarse and fine focusing drive, smoothness of coarse focusing drive being adjustable
- Mechanical stage 75×30 left/right with specimen holder, stage controls optionally on the right or left
- Space-saving, continuously adjustable illumination module optionally with halogen or LED source
- "Full-Köhler" or "Fixed-Köhler" Abbe condensers for brightfield, darkfield and phase contrast (for "Full-Köhler" design only)

- Backward inclined quadruple objective nosepiece with W 0.8 lens thread running on ball bearing
- Infinity-corrected "Plan-ACHROMAT" objectives with magnifications of 4x, 10x, 40x and 100x/Oil for brightfield, darkfield and phase contrast as well as for oil immersion applications (100x/Oil)
- Binocular tube or binocular phototube (50% vis, 50% doc) with ergonomically favorable tube angle of 30° with adjustable interpupillary distance and viewing height
- 10x adjustable eyepieces for field-of-view numbers 18 or 20, suitable for spectacle wearers

## 2.4 Objectives

The objectives are the optical heart of the microscope. The objectives may be labeled as follows:

Plan-ACHROMAT 10x/0.25 ∞/-.

Where:

10x Objective magnification, with a defined color ring on the objective being assigned to each magnification step (Carl Zeiss color code)

0.25 Numerical aperture

∞ Infinite mechanical tube length

– Usable with cover glass thickness  $D = 0$  or 0.17 mm

0 Usable without cover glass

or

0.17 Usable with cover glass thickness  $D = 0.17$  mm

Other labels:

Oil Oil immersion objective

Ph Phase contrast objective with green inscription

**Fig. 4 Objective**

Objective magnification multiplied by eyepiece magnification results in overall visual magnification, e.g.  $10 \times 10 = 100x$ .

Numerical aperture multiplied by 1000, e.g.  $0.25 \times 1000 = 250x$ , presents the maximum useful magnification; there is no resolution for further details above that limit.

When immersion objectives are used, the air between the cover glass and the objective is replaced by a liquid, which in most cases is immersion oil. The plastic oiler containing 5 ml of immersion oil is particularly suitable for this purpose.

Because of their short working distances, the 20x, 20x/Ph 2, 40x, 40x/Ph 2, 100x, 100x/Oil, and 100x/Ph 3 Oil objectives are equipped with resilient mounts for specimen protection.

### 3 START-UP AND OPERATION

#### 3.1 Starting up the microscope

##### 3.1.1 Setting up the microscope

Before installing and starting up the microscope, be sure to carefully read and observe the notes on instrument safety (see section 1).

Do not touch optical surfaces when unpacking the microscope to avoid fingerprints!


The microscope is supplied completely assembled and inclusive of accessories packed to commercial standards.

Additionally ordered components, such as sliders, transmitted-light module with illuminating mirror or reflected-light fluorescence illuminator, are delivered in separate packages and must still be mounted to the microscope.

- Remove the microscope from the transport case and place it onto the worktable.

Retain the original packaging for storage of the instrument in longer periods of non-use or for return to the manufacturer.


**Fig. 5 Starting up the microscope**


- Remove the plug-in power unit (Fig. 5/2) from its storage fixture on the back of the microscope stand.
- Replace the installed power outlet adapter by one of the supplied country-specific adapters (Fig. 5/4), if necessary. To this end, pull off the attached adapter and plug on the desired adapter.
- If you use a Primo Star iLED with reflected-light fluorescence illuminator, you will have to mount this module first, see section 3.1.2.
- If you use a battery supply unit, proceed according to section 3.1.3. If not, connect the plug-in power unit to a power outlet.
- If the plug-in power unit cannot be plugged into the chosen power outlet because of limited space, replace the power outlet adapter by the supplied IEC adapter (Fig. 5/3). This allows the plug-in power unit to be put flat onto the tabletop and connected to the power outlet through a country-specific appliance cable.

When using the IEC adapter, the plug-in power unit can be fixed for transportation purposes to the back of the microscope stand by means of the supplied two self-adhesive hooks and the Velcro® strip:


- Pull the Velcro® strip (Fig. 6/2) through the brackets of the hooks (Fig. 6/1).
- Stick the hooks (Fig. 6/1) together with the Velcro® strip to the recess (above on the right and left, respectively) at the back of the stand.
- Insert the plug-in power unit (Fig. 6/3) and close the Velcro® strip.

**Fig. 6** Fixing the plug-in power unit with IEC adapter to the stand

**Fig. 7** Mounting the magnetic field transmitter and the carrying handle

**Fig. 8** Placing the reflected-light illuminator onto the stand

### 3.1.2 Mounting the reflected-light fluorescence illuminator

#### Mounting the reflected-light illuminator onto the stand

- Remove the tube, see section 3.5.1. If the clamping screw is used to fix the tube to the stand, replace it by the Allen set screw.
- Insert the magnetic field transmitter (Fig. 7/1) for interlock as far as it will go into the stand at the rear strut and screw it down.
- Push the carrying handle (Fig. 7/3) from the back into the stand and fix it with two screws (Fig. 7/2).
- Insert the reflected-light illuminator (Fig. 8/1), in right angle and slightly inclined with its dovetail mount into the stand (Fig. 8/4).
- Position the reflected-light illuminator (Fig. 8/1) horizontally and turn it backwards with the dovetail mount in the stand, align it to the outer edges of the stand and tighten the set screw (Fig. 8/3).

If the outer edges of the reflected-light illuminator are not exactly aligned to the stand, it may be that the reflected-light illuminator cannot be switched on due to the interlock feature causing the interruption of the circuit to protect against emergent LED light.

- Place the tube onto the stand and tighten the clamping screw (Fig. 8/2), see section 3.5.1.

- Loosen the set screw (Fig. 9/6) on the back of the device, unplug the connector (Fig. 9/2) of the plug-in power unit (Fig. 9/5) and insert it into the connection socket (Fig. 9/1) of the reflected-light illuminator.
- Plug the connector (Fig. 9/3) of the reflected-light illuminator into the connection socket (Fig. 9/4) at the stand and tighten the set screw (Fig. 9/6).
- Fasten the cable of the reflected-light illuminator by pushing it into the cable clip (Fig. 9/7).
- Insert the plug-in power unit into a power outlet.


Fig. 9 Connecting the power supply unit

### Mounting the slider with yellow filter to the luminous-field diaphragm

- Unscrew the cover cap (Fig. 11/2) from the luminous-field diaphragm (Fig. 11/3).
- Turn the slider (Fig. 10/3) upside down. Loosen the clamping screw (Fig. 10/4) of the supporting ring and remove the supporting ring (Fig. 10/1) upward.
- Put the cover cap (Fig. 10/2), with its upper side facing downward, into the slider (Fig. 10/3).
- Insert the supporting ring (Fig. 10/1) into the slider and fix it with the clamping screw (Fig. 10/4).


Fig. 10 Inserting the cover cap into the slider


Fig. 11 Mounting the slider onto the luminous-field diaphragm


Fig. 12 Removing the yellow filter from the slider

- Turn the slider (Fig. 11/1) to its mounting position, put it on the luminous-field diaphragm (Fig. 11/3) and fasten it slightly, for the time being, by means of the cover cap (Fig. 11/2).
- Turn the slider (Fig. 11/1) to the desired position (for right-hand or left-hand operation) and fix it in this position by tightening the cover cap.

#### Removing the yellow filter from the slider

If you find the yellow filter (color conversion filter from blue (LED) to yellow) for transmitted-light observations bothersome, you may remove it as follows:

- Turn the cover cap of the luminous-field diaphragm (Fig. 11/2) to loosen it and remove it together with the slider (Fig. 11/1).
- Loosen the clamping screw (Fig. 12/6) and remove the filter holder (Fig. 12/5) downwards from the slider (Fig. 12/7).
- Remove the knurled ring (Fig. 12/4) with the color filter from the filter holder.
- Remove both O-rings (Fig. 12/1 and 2) from the knurled ring, push the yellow filter (Fig. 12/3) carefully out and keep it for future use.
- Put the knurled ring into the filter holder, insert both pieces together from below into the slider and fasten them with the clamping screw.

#### Using special eyecups with light protection

Our special eyecups with light protection (Fig. 17/1) are used for fluorescence applications if no darkroom is available. However, they **cannot be folded over** and are, therefore, not suitable for spectacle wearers.

For this reason, spectacle wearers should use the standard eyecups or fold-over eyecups.

### 3.1.3 Connecting the battery supply unit

#### Inserting batteries into the battery supply unit or remove them as follows:

- Loosen the four slotted screws of the battery supply unit.
- Remove the lid upward.
- Insert five commercial mono-cell batteries (D), type NiCd or NiMH, 1.2 V, with a capacity of 5000 mAh (min.) to 9000 mAh (max.), paying attention to correct polarity (see markings in the battery compartments).

Batteries of different types or different capacities must not be used together. Insert **rechargeable** batteries only.

- Push the changeover switch for the battery type (Fig. 3/1) to the correct position:  
**ON = NiMH; OFF = NiCd**
- Put the lid on again.
- Fasten it with the four slotted screws.

#### Primo Star iLED with reflected-light fluorescence illuminator:

- Pull the connector (Fig. 13/5) of the plug-in power unit (Fig. 13/6) out of the connection socket (Fig. 13/1) of the reflected-light illuminator and insert it into the connection socket of the battery supply unit (Fig. 13/4). The connector of the reflected-light illuminator (Fig. 13/3) is already connected to the connection socket of the microscope (Fig. 13/7).
- Insert the connector (Fig. 13/2) of the battery supply unit into the connection socket (Fig. 13/1) of the reflected-light illuminator.
- Insert the plug-in power unit (Fig. 13/6) into a power outlet.


Fig. 13 Connecting the battery supply unit

**Primo Star:**

- Loosen the set screw (Fig. 13/8) at the rear of the device and pull the connector (Fig. 13/5) of the plug-in power unit (Fig. 13/6) out of the connection socket (Fig. 13/7) of the microscope.
- Insert the connector of the plug-in power unit into the connection socket of the battery supply unit (Fig. 13/4).
- Insert the connector of the battery supply unit (Fig. 13/2) into the connection socket (Fig. 13/7) of the stand and tighten the set screw (Fig. 13/8).
- Insert the plug-in power unit (Fig. 13/6) into a power outlet.

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- Switch on the microscope with the rotary knob (Fig. 5/5) and adjust the illuminator to the desired intensity.

The selected intensity is indicated in five steps by the blue light-emitting diodes (Fig. 16/6) arranged on both sides of the stand.

- After finishing work, switch off the microscope with the rotary knob and cover it with the dust cover.
- The smoothness of the coarse focusing drive (Fig. 16/25 and Fig. 17/30) is factory-adjusted, but you can readjust it when required.

**Primo Star iLED:**

For transmitted-light applications:

Always turn the transmitted light / reflected light changeover switch first upward and then to the desired position. By using force to turn it downward the reflected-light illuminator will be damaged.

- Turn the transmitted light / reflected light changeover switch (Fig. 17/5) upward to the transmitted-light position (**Brightfield**).
- Switch on the transmitted-light illuminator using the rotary knob (Fig. 17/10) and adjust the desired illumination intensity.

The selected illumination intensity for transmitted light is indicated in five steps by the blue light-emitting diodes (Fig. 17/9) arranged on both sides of the stand.

For reflected-light applications (fluorescence):

Always turn the transmitted light / reflected light changeover switch first upward and then to the desired position. By using force to turn it downward the reflected-light illuminator will be damaged.

- Turn the transmitted light / reflected light changeover switch (Fig. 17/5) upward to the reflected-light position (**Fluorescence**).
- Switch on the reflected-light illuminator using the rotary knob (Fig. 17/6) and adjust the desired illumination intensity.

The pilot lamp (Fig. 17/25) at the front of the reflected-light illuminator lights with the reflected-light illuminator switched on. The brightness of the pilot lamp corresponds to the illumination intensity adjusted.

 The Primo Star iLED microscope with reflected-light fluorescence illuminator is provided with an interlock function, which causes the built-in LED of the reflected-light illuminator to be switched off as soon as the reflected-light illuminator is turned with regard to the stand or detached.

### Operation with battery supply unit:

When connecting the battery supply unit to a power outlet via the plug-in power unit of the microscope, the green power-on light **Ready** (Fig. 14/2) is on. The batteries are charged automatically.

The yellow charge indicator lamp **Chrg** (Fig. 14/3) lights as long as the charging process continues and goes out as soon as charging has been finished.

While charging is going on, the microscope can be used without any restriction. It is provided with voltage through the supply line.

If power supply to the supply line is interrupted or in case of power failure, the battery supply unit will automatically switch over to battery power. The power-on light **Ready** goes out.

Depending on the capacity of the batteries used, the microscope can be operated about 6 to 8 hours with battery power.

If the battery charge status has reached its critical value, the battery supply unit will switch off power supply to the microscope. In order to go on working and to charge the batteries, the unit must be connected to the supply line. Shortly before automatic disconnection of the rechargeable battery package (exhaustive discharge protection) the light begins to flash - now, at the latest, the battery package should be recharged at the supply line.

- The battery supply unit is switched on by pressing the **PowerOn** button (Fig. 14/1).
- After that, the reflected-light or the transmitted-light illuminator can be switched on at the microscope.

The battery supply unit need not be switched off. The unit will switch off automatically as soon as the reflected-light and the transmitted-light illuminators are switched off at the microscope.

- If the fuse has to be replaced, disconnect the plug-in power unit from the supply line and from the battery supply unit or pull off the connector of the plug-in power unit from the microscope.
- Unscrew the fuse holder (Fig. 15/1) from the battery supply unit and replace the defective fuse T4.0 A/H.
- Screw the fuse holder in and re-establish all cable connections again.


Fig. 14 Battery supply unit



Fig. 15 Battery supply unit, rear side

3.2 Controls of Primo Star (Full-Köhler and Fixed-Köhler)

English

A



Fig. 16 Controls of Primo Star

**Legend to Fig. 16:**

- 1 Eyepieces
- 2 Binocular body of the tube
- 3 Tube
- 4 Carrying handle
- 5 Plug-in power unit
- 6 Illumination-intensity indicators
- 7 Rotary knob for switching ON/OFF and adjusting the illumination intensity
- 8 Fine focusing drive (right side)
- 9 Coarse focusing drive (right side)
- 10 Control knob for X travel of mechanical stage
- 11 Control knob for Y travel of mechanical stage
- 12 Clamping screw for condenser
- 13 Transmitted-light illuminator, LED or HAL
- 14a Knurled ring for adjusting the luminous-field diaphragm (with Full-Köhler equipment only)
- 14b Luminous-field diaphragm (fixed in Fixed-Köhler equipment)
- 15a Centering screws for condenser on condenser carrier (in Full-Köhler equipment: knurled screws)
- 15b Centering screws for condenser on condenser carrier (in Fixed-Köhler equipment: Allen screws)
- 16a Abbe condenser, Full-Köhler
- 16b Abbe condenser, Fixed-Köhler
- 17 Objective
- 18 Microscope stage
- 19 Spring lever of specimen holder
- 20 Knurled ring of objective nosepiece
- 21 Lever for adjusting the aperture diaphragm of the condenser
- 22 Knurled knob for vertical adjustment of condenser
- 23 Coarse focusing drive (left side)
- 24 Fine focusing drive (left side)
- 25 Knurled ring for adjusting the smoothness of the coarse focusing drive

**Legend to Fig. 17:**

- 1 Special eyecups with light protection
- 2 Eyepieces
- 3 Binocular body of the tube
- 4 Tube
- 5 Transmitted light / reflected light changeover switch (Brightfield / Fluorescence)
- 6 Rotary knob for switching ON/OFF and adjusting the illumination intensity for reflected light
- 7 Carrying handle
- 8 Plug-in power unit
- 9 Illumination-intensity indicators for transmitted light
- 10 Rotary knob for switching ON/OFF and adjusting the illumination intensity for transmitted light
- 11 Fine focusing drive (right side)
- 12 Coarse focusing drive (right side)
- 13 Control knob for X travel of mechanical stage
- 14 Control knob for Y travel of mechanical stage
- 15 Clamping screw for condenser
- 16 Transmitted-light illuminator LED
- 17 Slider with yellow filter (with filter position for adapting the color temperature in transmitted light and with position for blocking the transmitted-light path in case of reflected-light fluorescence applications)
- 18 Luminous-field diaphragm (fixed)
- 19 Centering screws for condenser on condenser carrier
- 20 Abbe condenser, Fixed-Köhler
- 21 Objective
- 22 Microscope stage
- 23 Spring lever of specimen holder
- 24 Knurled ring of objective nosepiece
- 25 Pilot lamp for reflected-light fluorescence illuminator: lighting blue when switched on; brightness corresponds to intensity
- 26 Lever for adjusting the aperture diaphragm of the condenser
- 27 Knurled knob for vertical adjustment of condenser
- 28 Coarse focusing drive (left side)
- 29 Fine focusing drive (left side)
- 30 Knurled ring for adjusting the smoothness of the coarse focusing drive

### 3.3 Controls of Primo Star iLED (Fixed-Köhler) with reflected-light fluorescence illuminator


Fig. 17 Controls of Primo Star iLED

### 3.4 Operating the microscope

#### 3.4.1 Setting interpupillary distance and viewing height

- Rotate the eyepiece tubes symmetrically toward or away from one another to adjust the distance between the tubes to your individual interpupillary distance (Fig. 18).

The adjustment of the interpupillary distance is correct when you see only **one** round image while looking through the two eyepieces!

- Swivel the eyepiece tubes upward (Fig. 19/A) or downward (Fig. 19/B) to adjust the viewing height to your individual requirements.

#### 3.4.2 Compensating for ametropia and inserting the eyepiece pointer or eyepiece micrometer

The eyepieces (Fig. 20/3) are equipped with fold-over rubber eyecups (Fig. 20/1: pulled out; Fig. 20/2: folded over).

Both eyepieces are suitable for spectacle wearers. Additionally, they contain a focusing ring for the compensation of defective vision. The provided diopter scale serves to facilitate finding the correct setting.

When using the Primo Star iLED for fluorescence applications, the special eyecups with light protection can be used. However, they cannot be folded over and are not suitable for spectacle wearers.

If required, an eyepiece pointer or an eyepiece micrometer can be inserted in one eyepiece.

To this end, follow this procedure:

- Use Allen key SW 1 mm to loosen the set screw (Fig. 20/6) on the binocular body from below; remove the eyepiece.
- Unscrew the stop (Fig. 20/5) by hand from the eyepiece.
- Insert the eyepiece pointer (Fig. 20/4a) or the eyepiece micrometer (Fig. 20/4b) into the eyepiece (with the coated side facing your eyes). Screw in the eyepiece stop again.


Fig. 18 Setting the interpupillary distance


Fig. 19 Adjusting the viewing height



Fig. 20 Inserting the eyepiece pointer or eyepiece micrometer

- Put the eyepiece into the tube and fix it with the set screw.
- Turn the focusing ring of the eyepiece (Fig. 20/3) to focus on the wedge-shaped figure of the eyepiece pointer.
- Put the specimen onto the mechanical stage. Look at the specimen through the eyepiece with the eyepiece pointer and bring the microscopic image into focus using the focusing drive.
- When in the above-mentioned eyepiece both microscopic image and eyepiece pointer appear sharply defined, focus the image for the second eye by turning the focusing ring of the second eyepiece.

Having done so, both microscopic images inclusive of the eyepiece pointer are focused.

Afterwards, you should focus onto the specimen using the focusing drive only.

**Fig. 21 Adjusting transmitted-light brightfield**

- Look through one eyepiece of the binocular tube (Fig. 21/7) and bring the specimen into focus using the focusing drive (Fig. 21/2).
- Then, readjust image sharpness for the other eye, if necessary, by turning the eye lens of the adjustable eyepiece.
- Close the luminous-field diaphragm (Fig. 21/3) until it becomes visible (even if not in focus) in the field of view (Fig. 21/A).

### 3.4.3 Adjusting transmitted-light brightfield on the Full-Köhler microscope

- First, place a high-contrast specimen slide with the 0.17 mm cover glass being on top in the specimen holder of the mechanical stage. Fix the slide by means of the spring lever (Fig. 16/19).
- If the microscope stand is equipped with a phase or dark-field slider, pull this slider out to the left as far as it will go (phase slider with two phase contrast positions in mid-position).
- Adjust illumination intensity using the rotary knob (Fig. 21/1) on the microscope stand.
- Move the Abbe condenser up to the limit stop using the knurled knob (Fig. 16/22); set the control lever of the aperture diaphragm (Fig. 16/21) to mid-position.

The knurled knob for vertical adjustment of the condenser is on the left of the microscope if you use the mechanical stage 75x30 with drive on the right. It is on the right, if the stage controls are on the left.

- Rotate the 10x objective into the light path using the knurled ring (Fig. 21/6) of the nosepiece.

- Turn the knurled knob for vertical adjustment of the condenser (Fig. 16/22) until the edge of the luminous-field diaphragm appears sufficiently sharp (Fig. 21/B).
- Center the image of the luminous-field diaphragm using both centering screws (Fig. 21/5) of the condenser (Fig. 21/C). Then, open the diaphragm until it just disappears from the field of view (Fig. 21/D).
- To adjust the aperture diaphragm (contrast), remove one eyepiece from the tube and look through the tube with your naked eye. Swing lever (Fig. 21/4) to adjust the aperture diaphragm to approximately 2/3 ... 4/5 of the diameter of the exit pupil of the objective (Fig. 21/E). In most applications, this aperture diaphragm setting provides optimum contrast at almost ideal resolution, and is therefore the best compromise for the human eye.
- Insert the eyepiece again into the tube.

The size of the field of view and the objective aperture change after every objective change. Therefore, repeat the adjustment of luminous-field diaphragm and aperture diaphragm to obtain optimum results.

#### 3.4.4 Adjusting transmitted-light brightfield on the Fixed-Köhler microscope

The Primo Star Fixed-Köhler microscope is supplied factory-adjusted. Operation is restricted to a few manipulations.

- Place the specimen in the specimen holder of the mechanical stage.
- If the microscope stand is equipped with a dark-field slider, pull this slider to the left up to the lock-in position.
- When using the Primo Star iLED, turn the transmitted light / reflected light changeover switch to the transmitted light position (**Brightfield**) (turn it first fully upward). Push the slider containing the yellow filter with its filter position into the light path.
- Adjust the desired magnification by rotating the corresponding objective into the light path.
- Set the control lever of the condenser aperture diaphragm to the value of the selected magnification (**10x, 40x or 100x**).
- Focus onto the specimen using the focusing drive.
- Turn the rotary knob on the microscope stand to adjust the illumination intensity to a level pleasant for observation.

If the condenser had been removed (e.g. for installing the illuminating mirror), make sure to reinstall and center it by means of the two centering screws (for this, refer to section 3.5.5).


Fig. 22 Inserting the slider


Fig. 23 Centering the phase stop

### 3.4.5 Adjusting transmitted-light phase contrast or transmitted-light darkfield

The phase contrast method on the Primo Star microscope requires a stand with Full-Köhler equipment.

- First, adjust the microscope as you do for brightfield.
- Turn the nosepiece to position the phase-contrast objective (for Ph 1, Ph 2 or Ph 3) into the light path.
- Open the luminous-field diaphragm (Fig. 22/3) on the stand and open the aperture diaphragm by means of the control lever (Fig. 22/1) on the Abbe condenser (Fig. 22/2).
- Insert the slider fitting the objective used, with one (Fig. 22/6) or two (Fig. 22/9) phase contrast positions (Ph).
  - Slider Ph 1, Ph 2 or Ph 3 (Fig. 22/6): Remove the screw (Fig. 22/7). Push the slider from the left into the Abbe condenser and screw in the screw (Fig. 22/7) again.
  - Slider Ph 1 / H / Ph 2 (Fig. 22/9): Unscrew the handle (Fig. 22/8) at the right side. Push the slider (Fig. 22/9) from the left into the Abbe condenser (the inscriptions Ph 1 and Ph 2 are readable upright and laterally correct from the front). Screw in the handle again.


#### Caution: Risk of pinched fingers

If the slider Ph 1 / H / Ph 2 is used, your fingers might be pinched between the slider and the control knobs of the microscope stage. This risk can be prevented by removing the slotted screw (on the underside of the condenser) from the middle bore hole and screwing it into the bore hole to the left or right of it (depending on whether the stage is operated from the right or from the left).

- First, loosen the clamping screw (Fig. 16/12) on the condenser carrier and pull the condenser out to the front (if necessary, lower the condenser carrier using the knurled knob (Fig. 16/22)).
- Unscrew the slotted screw (Fig. 24) and screw it into the corresponding bore hole to the right or left of it. Insert the condenser again.
- Push the slider to the right until perceiving the lock-in position (to the right or left in the case of the slider with two phase contrast positions) to move the phase stop into the light path.

 The filter ( $\varnothing$  22 mm) in the mid-position of the slider Ph 1 / H / Ph 2 has been factory-installed. After removing the supporting ring it can be exchanged or removed.

**Fig. 24** Slotted screw on the underside of the condenser

- Open the aperture diaphragm of the condenser completely using its control lever (left limit stop).
- Adjust the illumination intensity as needed.
- Check the centering of the phase stop according to the drawing in Fig. 23. To this end, remove one eyepiece and replace it by the diopter.
- Center the phase stop, if necessary, (Fig. 23/A) by turning the two adjusting screws of the slider (Fig. 22/5) by means of the two Allen keys SW 1.5 (Fig. 22/4) until the phase stop image corresponds to that shown in (Fig. 23/B).
- Afterwards, replace the diopter by the eyepiece again.


 For darkfield applications, use the darkfield slider in place of the phase-contrast slider.

### 3.4.6 Adjusting reflected-light fluorescence

To carry out examinations according to the reflected-light fluorescence method, the Primo Star iLED microscope stand with reflected-light fluorescence illuminator is required.

The Primo Star iLED microscope with reflected-light fluorescence illuminator is provided with an interlock function, which causes the built-in LED of the reflected-light illuminator to be switched off as soon as the reflected-light illuminator is turned with regard to the stand or detached.

**Fig. 25** Adjusting reflected-light fluorescence

- First, adjust the microscope as you do for brightfield, see section 3.4.4.
- Turn the nosepiece to position the objective for fluorescence application into the light path (e.g. objective 40x).
- Set the transmitted light / reflected light changeover switch (Fig. 25/1) to the reflected light position (**Fluorescence**) (first turn it fully upward).
- Switch on the reflected-light LED using the rotary knob (Fig. 25/2) of the reflected-light illuminator and adjust the illumination intensity to a level pleasant for observation. The pilot lamp (Fig. 25/4) at the front of the reflected-light illuminator shows blue. The brightness of the pilot lamp corresponds to the illumination intensity adjusted for reflected light.
- Focus on the specimen using the focusing drive.
- To avoid interfering fluorescences (caused by the transmitted-light LED), push the blocking position of the slider (Fig. 25/3) into the light path.

If the slider is not available, about 90% of the interfering fluorescences can be eliminated by closing the condenser aperture diaphragm.

Our special eyecups with light protection (Fig. 25/5) can be used for microscopic examinations in rooms which are not darkened. They are, however, not suitable for spectacle wearers and must not be folded over because, otherwise, the required dimensional stability would get lost.

### 3.5 Converting the microscope

Unplug the plug-in power unit from line power before converting the microscope.

#### 3.5.1 Changing the tube

- Loosen the clamping screw (Fig. 26/2), turn the installed tube (Fig. 26/1) clockwise by approximately 90° (Fig. 26/A) and detach it on the right side upward (Fig. 26/B).

For reasons of space, the tube may also be clamped by the Allen set screw supplied with the tube.

- Insert the tube to be installed, with the eyepieces pointing to the right and its dovetail mount slightly inclined, under the two holding elements (Fig. 26/4) into the stand.
- Then, place the tube horizontally onto the stand. In doing so, the groove on the underside of the tube must be located above the third holding element (Fig. 26/3) of the stand.
- Turn the tube counterclockwise by 45° (eyepieces pointing frontward), align it to the stand and tighten the clamping screw (Fig. 26/2).

For space-saving storage of the microscope (e.g. in a cabinet) the tube may also be rotated by 180° to the back.


Fig. 26 Changing the tube




Fig. 27 Inserting a color filter

#### 3.5.2 Inserting color filters

- Move the condenser carrier fully up by turning the knurled knob for vertical adjustment (Fig. 16/22).
- Unscrew the cover cap (Fig. 27/3) from the luminous-field diaphragm (Fig. 27/1).
- Put the desired filter – yellow, green or blue – (Fig. 27/2) onto the mounting surface of the luminous-field diaphragm and screw the cover cap back into place.

**Fig. 28** Replacing the 6V 30W halogen lamp or the LED illumination module

**Fig. 29** Changing an objective

### 3.5.3 Replacing the 6V 30W halogen lamp or the LED illumination

Unplug the plug-in power unit from the power outlet and allow for a sufficient cool-down time of the 6V 30W halogen lamp before you replace it.

- Loosen both fastening screws (Fig. 28/3) of the illumination module (Fig. 28/1 or 2). In doing so, slightly press the screws against the spring and turn them by 90°: Turn left screw clockwise and right screw counterclockwise.
- Pull the illumination module (Fig. 28/1 or 2) out of the stand.
- If the stand is equipped with a 6V 30W halogen lamp (Fig. 28/2), remove the halogen lamp (Fig. 28/4) from the lamp carrier and insert the new halogen lamp. Do not touch the new lamp with naked fingers as this will reduce lamp life.
- If the stand is equipped with an LED source, replace the complete illumination module including the LED (Fig. 28/1).
- Push the illumination module (Fig. 28/1 or 2) back into the stand and fasten it by means of the two screws (Fig. 28/3). In doing so, slightly press the screws against the spring and turn them by 90°: Turn left screw counterclockwise and right screw clockwise.

### 3.5.4 Changing objectives

- Turn the focusing drive to move the mechanical stage fully down.
- Turn the nosepiece (Fig. 29/3) to move the objective to be changed (Fig. 29/1) into a lateral position.
- Unscrew the objective using the rubber strip included in the delivery and remove it downward.
- Screw the desired objective (Fig. 29/2) hand-tight into the nosepiece as far as it will go.
- If you intend to insert the objective in a lens mount not used so far, remove the dust cover from the corresponding mount of the nosepiece.

### 3.5.5 Installing/removing the mirror

The mirror serves to illuminate the specimen if no power outlet is available.

It can be used only in combination with the Primo Star Fixed-Köhler Microscope. For this purpose, it is necessary to remove the condenser and the condenser carrier insert.

#### Installing the mirror:

- Remove the cover cap (Fig. 30/6) from the luminous-field diaphragm (Fig. 30/5).
- Unscrew the stop screw (Fig. 30/8) and lower the condenser carrier fully down using the knurled screw for vertical adjustment.
- Using an Allen key, loosen the clamping screw (Fig. 30/2) of the condenser and the centering screws (Fig. 30/3) of the condenser carrier (Fig. 30/4) so that the condenser (Fig. 30/7) can be removed frontward. Remove the condenser (Fig. 30/7).
- Unscrew the clamping screw (Fig. 30/2 or Fig. 31/1) from the condenser carrier insert (Fig. 30/1 or Fig. 31/3).
- Push the condenser carrier insert backward against the spring and inclining it remove it upward from the condenser carrier (Fig. 30/4).
- Insert the mirror (Fig. 31/5) from top through the opening of the condenser carrier and place it onto the mount of the luminous-field diaphragm (Fig. 31/2). Take care that the mirror rests horizontally on the mount.
- Rotate and incline the mirror until the daylight is reflected homogeneously into the light path.

#### Removing the mirror:

Removal of the mirror is done in reversed order.


Fig. 30 Removing/installing the condenser


Fig. 31 Installing/removing the mirror


Fig. 32 Installing a camera

### 3.5.6 Installing a camera

You can install a digital camera, a video camera or a compact digital camera to your choice on the microscope via the available camera adapters (see section 2.1).

- First, install the phototube (Fig. 32/7), if required (see section 3.5.1).
- Loosen the clamping screw (Fig. 32/8) and remove the dust cap from the phototube.

#### Installing a compact digital camera

- Sliding mount (Fig. 32/4), thread adapter ring M37/52 (Fig. 32/3) and lens mount (Fig. 32/6) are supplied pre-assembled as **Digital Camera Adapter P95 M37/52x0.75**. The opposite drawing shows this unit dismantled. Furthermore, you may unscrew the M37/52 thread adapter (Fig. 32/3) from the sliding mount (Fig. 32/4) (not shown) so that you can also mount cameras with M37 thread.
- Mount the adapter ring (Fig. 32/2) (supplied with the camera) to the camera (Fig. 32/1) (see the operating instructions of the camera).
- Screw the unit consisting of sliding mount (Fig. 32/4), thread adapter ring M37/52 (Fig. 32/3) and lens mount (Fig. 32/6) into the adapter ring (Fig. 32/2).
- Insert the camera with adapter into the phototube as far as it will go. Align it and fasten it with the clamping screw (Fig. 32/8).
- Depending on the microscope equipment or the camera used, it may be necessary to optimize the distance between camera lens and lens mount (Fig. 32/6) (see double-headed arrow). This will be necessary in particular, if it is impossible to obtain an unvignetted image in any of the zoom positions of the camera lens. To this end, adjust the camera as follows:
  - Switch off autofocus.
  - Set the object distance to  $\infty$ .
  - Set aperture-priority auto exposure mode.
  - Choose an aperture as large as possible (i.e. small aperture value!).
 Not all cameras provide these options. Please consult the Operating Manual of the camera used.
- Loosen the set screw (Fig. 32/5).
- Vary the camera lens/lens mount distance gradually, i.e. displace the sliding mount with camera on lens mount in defined steps.
- Zoom the camera lens through from wide angle (W) to tele position (T).

- Carry out this test until the image is format filling without masking or vignetting.
- Retighten the set screw (Fig. 32/5).

With camera/adapter combinations that have not been recommended expressly by Zeiss, it may be quite impossible to obtain an unvignetted image.

### Installing a digital camera or video camera with C-mount thread

Cameras with C-mount thread are to be connected to the phototube of the microscope by means of the camera adapter P95-C 2/3" 0.65x or P95-C 1/2" 0.5x (Fig. 32/9).

- Put the camera (Fig. 32/10) together with the corresponding camera adapter (Fig. 32/9) into the phototube as far as it will go, align it and fasten it by means of the clamping screw (Fig. 32/8).

## 4 CARE AND TROUBLESHOOTING

### 4.1 Instrument care

Care of the microscope is restricted to the following operations:

- Cover the instrument with the dust cover after every use.
- Do not install the microscope in a humid room; maximum humidity < 75%.
- Cover open tubes with the dust caps.
- Remove dust and loose dirt from visible optical surfaces with a brush, blower brush, cotton bud, optics cleaning tissue, or a cotton cloth.
- Remove water-soluble dirt (coffee, cola, etc.) by blowing on it and subsequent wiping it off with a dust-free cotton cloth or a cloth moistened with water to which you may also add a mild detergent.
- Wipe off stubborn oily or fatty dirt (immersion oil, fingerprints) with a cotton bud or a dust-free cotton cloth and the optics cleaning solution L.

The cleaning solution consists of 90 vol% gasoline and 10 vol% isopropanol (IPA). Its individual constituents are also known as:

Gasoline:	Surgical spirit, petroleum ether
Isopropanol:	2-propanol, dimethyl carbinol, 2-hydroxypropane

Clean optical surfaces by polishing in circles starting in the middle and moving to the edges using slight pressure only.

Disconnect the plug-in power unit from the supply line before cleaning it. Make sure that no moisture gets into the plug-in power unit.

For the use in warm and humid climatic zones, all optical components of the microscope are provided with protection against fungus attack.

## 4.2 Troubleshooting

Problem	Cause	Remedy
The field of view is not completely visible.	Nosepiece with objective has not been switched into click-stop position.	Switch nosepiece with objective into click-stop position.
	Condenser has not been set correctly.	Set condenser correctly.
	The aperture diaphragm has not been adjusted correctly.	Adjust aperture diaphragm correctly.
	The luminous-field diaphragm has not been adjusted correctly.	Adjust luminous-field diaphragm correctly.
	The filter has not been inserted correctly into the filter mount.	Insert filter correctly into filter mount.
Low resolving power, poor image contrast	The aperture diaphragm has not been opened to correct size.	Set opening of aperture diaphragm to correct size.
	Condenser not focused correctly.	Focus the condenser.
	Wrong cover glass thickness selected for use of transmitted-light objectives corrected for 0.17mm cover glass.	Use standard 0.17 mm cover glass.
	Use of no or non-specified immersion oil with immersion objective	Use the supplied immersion oil.
	Air bubbles in immersion oil	Remove the bubbles by applying new oil or moving the objective to and fro.
	Immersion oil on the front lens of a dry objective Dirt or dust on the optical surfaces of objectives, eyepieces, condenser, or filters	Clean the front lens of the dry objective. Clean the respective optical components.
Greater focus differences after changing the objective	The adjustable eyepieces have not been set correctly.	Set the adjustable eyepieces to the corresponding defective vision.
The 6V 30W halogen lamp or the LED source does not light up although the microscope has been switched on.	Power plug has not been plugged into power outlet.	Connect power plug to power outlet.
	6V 30W halogen lamp or LED source is defective.	Replace the defective 6V 30W halogen lamp or the LED module.
The 6V 30W halogen lamp is flickering, its light intensity unstable, the illumination inhomogeneous	End of average service life of 6V 30W halogen lamp has been reached.	Replace the 6V 30W halogen lamp.
	Power cable incorrectly installed or broken.	Connect the power cable correctly or replace it.
	The pins of the 6V 30W halogen lamp have not been inserted correctly into the lamp holder.	Insert the pins of the 6V 30W halogen lamp correctly into the lamp holder.
	The pins of the 6V 30W halogen lamp have not been inserted symmetrically into the lamp holder.	Insert the pins of the 6V 30W halogen lamp symmetrically into the lamp holder.
Stage comes down by itself, image focus is unstable	Adjusted torque of coarse focusing drive is too low.	Increase the torque of the coarse focusing drive to make motion stiffer.

Problem	Cause	Remedy
Transmitted-light illuminator of Primo Star iLED cannot be switched on.	Transmitted light / reflected light changeover switch is in reflected light ( <b>Fluorescence</b> ) position.	Turn transmitted light / reflected light changeover switch to transmitted light ( <b>Brightfield</b> ) position.
Reflected-light illuminator of Primo Star iLED cannot be switched on.	Transmitted light / reflected light changeover switch is in transmitted light ( <b>Brightfield</b> ) position.	Turn transmitted light / reflected light changeover switch to reflected light ( <b>Fluorescence</b> ) position.
	Reflected-light illuminator not exactly enough aligned to microscope stand or displaced, causing interruption of power supply by interlock function.	Align reflected-light illuminator exactly to microscope stand and tighten clamping screw sufficiently.

### 4.3 Changing the LED module in the reflected-light fluorescence illuminator

#### ATTENTION

The LED module in the reflected-light fluorescence illuminator has to be changed in any case by an authorized service technician.

The Primo Star iLED microscope with reflected-light fluorescence illuminator is provided with an interlock function to switch off the installed LED of the reflected-light illuminator as soon as the reflected-light illuminator will be turned with regard to the stand or detached or if the lid of the reflected-light illuminator will be removed.

- Disconnect the reflected-light illuminator from the power supply.
- Detach the tube (see also section 3.5.1). To this end, loosen the knurled screw, turn the tube counterclockwise by approximately 90° and lift it out of the dovetail mount.
- Loosen the three fastening screws (Fig. 33/1) of the lid of the reflected-light illuminator (2.5 mm Allan key).
- Disconnect the plug-in connector of the LED module connection cable (Fig. 34/3) on the circuit board. To this end, press the safety hook against the connector to release it and pull out the connector.
- Unscrew the two fastening screws of the LED holder (Fig. 34/1) completely. Remove holder (Fig. 34/2) and screws.
- Remove the LED module (Fig. 34/4).
- Insert the new LED module. Because of limited space, the flat side of the LED module must point to the circuit board.



Fig. 33 Lid of reflected-light illuminator

**Fig. 34 Reflected-light LED module**

- Attach the holder (Fig. 34/2) again, tighten the screws (Fig. 34/1) only slightly, still allowing displacement of the LED module within the guide.
- Shift the LED module frontward (operator side) against the stop and fix the screws (Fig. 33/1) handtight. This step is important to re-establish the optimally adjusted optical configuration.
- Plug the connector (Fig. 34/3) again into the circuit board.
- Put on the lid and tighten the three fastening screws (Fig. 33/1) after shifting the lid in direction of arrow to the front right corner (Fig. 33) to realign the correct optical axis.
- Attach the tube again (see also section 3.5.1).

## 5 APPENDIX

### 5.1 Technical data

#### Dimensions (width x depth x height)

Stand with binocular tube	approx. 190 x 410 x 395 mm
Stand with phototube	approx. 190 x 425 x 395 mm
With tube / phototube turned by 180°	approx. 190 x 375 x 395 mm
Stand with reflected-light fluorescence illuminator	approx. 190 x 410 x 449 mm

#### Weight

Primo Star with phototube	approx. 8.2 kg
Primo Star iLED with reflected-light fluorescence illuminator and phototube	approx. 9.6 kg

#### Ambient conditions

Transportation (in packaging): Permissible ambient temperature	-40 to +70 °C
Storage: Permissible ambient temperature	+10 to +40 °C
Permissible air humidity (no condensation)	max. 75 % at 35 °C
Operation: Permissible ambient temperature	+10 to +40 °C
Permissible air humidity (no condensation)	max. 75 % at 35 °C
Atmospheric pressure	800 hPa to 1060 hPa

#### Operating data

Protection class	II
Protection type	IP20
Electrical safety	in compliance with DIN EN 61010-1 (IEC 61010-1) including CSA and UL directives
Pollution degree	2
Overtoltage category	II
Radio interference suppression	in accordance with EN 61326
Line voltage	100 to 240 V (±10 %) wide-range input power supply, i.e. voltage setting of the instrument need not be changed!
Line frequency	50 / 60 Hz
Power consumption	70 VA; secondary voltage of external power supply 12 V
Plug-in power unit output	12 V DC; max. 2.5 A
Microscope 12 V / 6 V DC	adjustable from 1.5 V to 6 V
LED class of complete device	3B

#### Light sources

Halogen lamp	HAL 6 V, 30 W
Adjustability of light source	continuous, from 1.5 to 6 V DC
Color temperature at 6V	2800 K
Luminous flux	280 lm
Average service life	1000 h
Luminous area	1.5 x 3 mm

LED illumination	white light LED, peak wavelength 440 nm, LED class 2
Constant, brightness-independent color temperature of	7480 K
Homogeneous field illumination	20 mm diameter
Suitable for objectives with magnifications of	4x to 100x
Analogous brightness adjustment from	approx. 15 to 100 %
LED modules (reflected-light fluorescence illuminator)	max. 40 mW, 365 - 625 nm; LED class 3B

**Battery supply unit (accessory)**

Batteries	fuses according to IEC 127 T4.0 A/H
Type	mono-cell (D) - commercially available, NiCd or NiMH, 1.2 V
Capacity	minimum 5000 to max. 9000 mAh
Number per battery supply unit	5 batteries
Operational lifetime	several hours, depending on the capacity of the batteries

**Optical/mechanical data**

Stand with stage focusing	
With coarse focusing drive	45 mm/rev.
With fine focusing drive	0.5 mm/rev.
Total stage lift	15 mm
Objective change	manual via quadruple objective nosepiece
Objectives	infinity-corrected objective range with W 0.8 mounting thread
Eyepieces	30 mm tube size
With field-of-view number 18	PL 10x/18 Br. foc.
With field-of-view number 20	PL 10x/20 Br. foc.
Specimen stage	mechanical stage 75x30 right/left
Dimensions (width x depth)	140 x 135 mm
Stage travel (X x Y)	75 x 30 mm
Coaxial drive	optionally right or left
Vernier scales	readable from the right
Specimen holder	with spring lever, left
Abbe condenser 0.9/1.25; Fixed-Köhler	for $V_{obj}$ 4x to 100x
Abbe condenser 0.9/1.25; Full-Köhler	for $V_{obj}$ 4x to 100x
Binocular tube 30°/20	
Maximum field-of-view number	20
Interpupillary distance	adjustable from 48 to 75 mm
Tube angle	30°
Viewing height	380 to 415 mm
Viewing port	tube factor 1x
Binocular phototube 30°/20	
Maximum field-of-view number	20
Interpupillary distance	adjustable from 48 to 75 mm
Tube angle	30°
Viewing height	380 to 415 mm
Viewing port	tube factor 1x
Photo/video port	tube factor 1x, 60 mm mount
Invariable splitting ratio	50% vis / 50 % doc
Illuminating mirror	with plane surface and spherical surface with $f' = 75$ mm

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